

**Joint ICTP-IAEA Workshop on Monte Carlo Radiation Transport
and Associated Data Needs for Medical Applications**

28 October – 8 November 2024

ICTP, Trieste, Italy

Lecture 19

DOSXYZnrc calculations with CT input

Blake Walters

Metrology Research Centre
National Research Council Canada



Government
of Canada


Gouvernement
du Canada



DOSXYZnrc phantom definition: CT input

- DOSXYZnrc can use phantoms derived from CT data sets allowing simulation in realistic anthropomorphic phantoms.
- Voxelized phantoms from CT data are obtained using the stand-alone application `ctcreate`
- `ctcreate` supports CT data in the following formats: DICOM, ADAC Pinnacle, and CADPLAN.
- A tool for converting the AAPM CT format into Pinnacle format is also available

DOSXYZnrc phantom definition: GUI CT input



Input full name of
.egsphant file

- Currently zeroes dose in all voxels with $\rho \leq 0.44 \text{ g/cm}^3$
- This can be changed in `dosxyznrc.mortran`
- For display only


Recommend supressing to prevent
Mbyte -- Gbyte .egs1st file

CT scan to *.egsphant tool: **ctcreate**


- Stand alone MORTRAN code for converting CT image data to .egsphant files
- Resides in \$OMEGA_HOME/progs/ctcreate
- Can be used with DICOM and Pinnacle format CT data
- Run using: ctcreate [inputfile]

ctcreate: flow chart


A flowchart for use of CT data with **ctcreate** and **DOSXYZnrc**



ctcreate: Subset and voxel selection



DOSXYZnrc Phantom (.egsphant)*



ctcreate: Input file format

- full name of header file (Pinnacle) ← or DICOM
 - name of file listing image files (DICOM) → /Users/blkewalters/EGSnrc/egs_home/dosxyznrc/image_001.header
 - CT conversion ramp →
 - Input "0,0" on first line for default ramp
- ```
Pinnacle ← or DICOM
7.0, 25.0, 23.50, 46.0, -9.0, 5.0 ← X, Y, Z limits of subset of CT data
.5, .5, .5 ← X, Y, Z voxel dimensions for phantom
4, 1 ← no. of media, HUlower(1)
AIR700ICRU ← medname(1)
50, 0.001, 0.044, 1.0 ← HUupper(1), rholower(1), rhoupper(1), dummy(1)
LUNG700ICRU
300, 0.044, 0.302, 1.0
ICRUTISSUE700ICRU
1125, 0.302, 1.101, 1.0
ICRPBONE700ICRU
3000, 1.101, 2.088, 1.0
~
```

## ctcreate: Input DICOM CT data

```
DICOM
/Users/blakewalters/EGSnrc/egs_home/dosxyznrc/mouse_CT_list
-25.0, 25.0, -15.0, 35.0, 0.0, 3.0
0.025, 0.025, 0.1
5, -1000
Air
-850, 0.0012, 0.02, 1.0
Lung
-400, 0.02, 0.6, 1.0
softTissue1
0, 0.6, 1.0, 1.0
softTissue2
227, 1.0, 1.05, 1.0
Bone
3723, 1.05, 1.82, 1.0
```

**A file containing a list of DICOM  
image files in order of  
increasing Z (slice position)**


```
/Users/blakewalters/EGSnrc/egs_home/dosxyznrc/mouse_CT/slice_0000.dcm
/Users/blakewalters/EGSnrc/egs_home/dosxyznrc/mouse_CT/slice_0001.dcm
/Users/blakewalters/EGSnrc/egs_home/dosxyznrc/mouse_CT/slice_0002.dcm
/Users/blakewalters/EGSnrc/egs_home/dosxyznrc/mouse_CT/slice_0003.dcm
/Users/blakewalters/EGSnrc/egs_home/dosxyznrc/mouse_CT/slice_0004.dcm
/Users/blakewalters/EGSnrc/egs_home/dosxyznrc/mouse_CT/slice_0005.dcm
/Users/blakewalters/EGSnrc/egs_home/dosxyznrc/mouse_CT/slice_0006.dcm
/Users/blakewalters/EGSnrc/egs_home/dosxyznrc/mouse_CT/slice_0007.dcm
/Users/blakewalters/EGSnrc/egs_home/dosxyznrc/mouse_CT/slice_0008.dcm
/Users/blakewalters/EGSnrc/egs_home/dosxyznrc/mouse_CT/slice_0009.dcm
/Users/blakewalters/EGSnrc/egs_home/dosxyznrc/mouse_CT/slice_0010.dcm
/Users/blakewalters/EGSnrc/egs_home/dosxyznrc/mouse_CT/slice_0011.dcm
/Users/blakewalters/EGSnrc/egs_home/dosxyznrc/mouse_CT/slice_0012.dcm
/Users/blakewalters/EGSnrc/egs_home/dosxyznrc/mouse_CT/slice_0013.dcm
/Users/blakewalters/EGSnrc/egs_home/dosxyznrc/mouse_CT/slice_0014.dcm
/Users/blakewalters/EGSnrc/egs_home/dosxyznrc/mouse_CT/slice_0015.dcm
/Users/blakewalters/EGSnrc/egs_home/dosxyznrc/mouse_CT/slice_0016.dcm
/Users/blakewalters/EGSnrc/egs_home/dosxyznrc/mouse_CT/slice_0017.dcm
/Users/blakewalters/EGSnrc/egs_home/dosxyznrc/mouse_CT/slice_0018.dcm
/Users/blakewalters/EGSnrc/egs_home/dosxyznrc/mouse_CT/slice_0019.dcm
/Users/blakewalters/EGSnrc/egs_home/dosxyznrc/mouse_CT/slice_0020.dcm
/Users/blakewalters/EGSnrc/egs_home/dosxyznrc/mouse_CT/slice_0021.dcm
/Users/blakewalters/EGSnrc/egs_home/dosxyznrc/mouse_CT/slice_0022.dcm
/Users/blakewalters/EGSnrc/egs_home/dosxyznrc/mouse_CT/slice_0023.dcm
```



## ctcreate: DICOM CT data

- Uses the C subroutine, `ReadCT_DICOM.C`, which is linked with `ctcreate` during compilation
- Have to set macro `MAX_SLICES` large enough for number of image slices
- Gives you the option of applying an offset to all HU numbers – may be useful depending on definition of CT conversion ramp
- Is fairly general (give or take some minor tweaking) and may provide a good place to start if you do not have your own routine for converting DICOM images

## ctcreate: Default conversion ramp



# ctcreate: Output

```
Input the format of your CT data
:
Pinnacle

Input the full name of the header file for the CT data
: /Users/blakewalters/EGSnrc/egs_home/dosxyznrc/image_001.header
```

```

CT Phantom has been chosen and reading
headerfile information.
```

```
Header File Name -> /Users/blakewalters/EGSnrc/egs_home/dosxyznrc/image_001.header
X Pixel Number -> 256
Y Pixel Number -> 256
Z Pixel Number -> 28
X Pixel Size -> 0.123 (cm)
Y Pixel Size -> 0.123 (cm)
Z Pixel Size -> 1.000 (cm)
X Pixel Offset -> 0.000 (cm)
Y Pixel Offset -> 17.400 (cm)
Z Pixel Offset -> -9.000 (cm)
```

- The macros `$CTIMAX`, `$CTJMAX`, `$CTKMAX` in `ctcreate.mortran` must be large enough to accommodate these dimensions
- If not, change them and then recompile `ctcreate`

```
CT data goes from x = 0.00000 - 31.36000 cm
 y = 17.40000 - 48.76000 cm
 z = -9.00000 - 19.00000 cm
```

```

Reading in the CT data from the image file.
/Users/blakewalters/EGSnrc/egs_home/dosxyznrc/image_001.img
Done reading in the image file data.
```

```

Determining if data needs to be byte
swapped.
```

```
Byte order of file (0) is same as byte order of this machine 1234
No Byte swap necessary, continuing.
```

```
Warning: 3226 pixels had Pinnacle CT no. > the maximum value of 2000
supported by Pinnacle format (or indicated in header).
```

Useful output  
about CT  
data read in

## ctcreate: Output (cont.)

CT Volume subset selection.

Please enter the positions of limiting  
planes (cm):

```
xctsubmin,xctsubmax,yctsubmin,yctsubmax,zctsubmin,zctsubmax
: 7.0000 25.0000 23.5000 46.0000 -9.0000 5.0000
```

The voxel index limits are as follows:

I Limits -> i= 58 to i= 205

J Limits -> j= 50 to j= 234

K Limits -> k= 1 to k= 15

```
xctsubmin,xctsubmax,yctsubmin,yctsubmax,zctsubmin,zctsubmax (cm)
```

after adjustment to fit integer no. of voxels

```
: 6.9825 25.1125 23.4025 46.0651 -9.0000 6.0000
```

**Note: subset boundaries adjusted to  
fit an integer no. of CT voxels**

Resample CT data for dosxyznrc

Input the x,y,z dimensions (cm) of the dosxyznrc voxel on one line

```
(min= 0.14164 x 0.17705 x 0.11719 cm) ←
```

```
: 0.50000 0.50000 0.50000
```

**Min. phantom voxel dimensions determined by  
\$IMAX, \$JMAX, \$KMAX, in  
\$HEN\_HOUSE/user\_codes/dosxyznrc/dos  
xyznrc\_user\_macros.mortran**

New X voxel thickness -> 0.50

New Y voxel thickness -> 0.50

New Z voxel thickness -> 0.50

New number X voxels -> 36

New number Y voxels -> 45

New number Z voxels -> 30

Final x,y,z dimensions of dosxyznrc voxels in cm (adjusted so that an integer  
number fit exactly on the CT data): 0.50361 0.50361 0.50000 ←

Calculating bounds and new CT values

**Note: phantom voxel dimension  
adjusted to fit an integer no.  
into the CT subset chosen**

## ctcreate: Output (cont.)

The CT-Density Ramp

Number of media (max 7), min. CT number of ramp  
(0,0 if you want to use the hard-wired ramp function): 4 1

Medium 1 : AIR700ICRU  
CT no. upper bound, density lower bound (g/cm<sup>3</sup>),  
density upper bound (g/cm<sup>3</sup>)--all on one line  
: 50 0.00100 0.04400

Medium 2 : LUNG700ICRU  
CT no. upper bound, density lower bound (g/cm<sup>3</sup>),  
density upper bound (g/cm<sup>3</sup>)--all on one line  
: 300 0.04400 0.30200

Medium 3 : ICRUTISSUE700ICRU  
CT no. upper bound, density lower bound (g/cm<sup>3</sup>),  
density upper bound (g/cm<sup>3</sup>)--all on one line  
: 1125 0.30200 1.10100

Medium 4 : ICRPBONE700ICRU  
CT no. upper bound, density lower bound (g/cm<sup>3</sup>),  
density upper bound (g/cm<sup>3</sup>)--all on one line  
: 3000 1.10100 2.08800

Writing CT phantom data into image\_001.egsphant to be read by dosxyznrc.

← CT conversion ramp

← Appends .egsphant to image file  
name & writes it out to the current  
directory

# Phantom definition: **\*.egsphant** format


| block | data                                             | description                                                                                                                                                                                                              |
|-------|--------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1     | nmed                                             | no. of media in simulation                                                                                                                                                                                               |
| 2     | medname(1)<br>medname(2)<br>...<br>medname(nmed) | <ul style="list-style-type: none"> <li>List of medium names</li> <li>Any medium surrounding phantom (dsurround) must be included here</li> <li>correspond to media in PEGS file or defined in pegsless inputs</li> </ul> |
| 3     | smax(1) ... smax(nmed)                           | a single line of inputs no longer used                                                                                                                                                                                   |
| 4     | nx ny nz                                         | no. of voxels in X, Y, Z direction                                                                                                                                                                                       |
| 5     | (xbound(i),i=1,nx+1)                             | X voxel boundaries                                                                                                                                                                                                       |
| 6     | (ybound(j),j=1,ny+1)                             | Y voxel boundaries                                                                                                                                                                                                       |
| 7     | (zbound(k),k=1,nz+1)                             | Z voxel boundaries                                                                                                                                                                                                       |
| 8     | ((med(i),i=1,nx),j=1,ny),k=1,nz)                 | <ul style="list-style-type: none"> <li>Voxel medium indices in I1 (single integer)</li> <li>&lt;return&gt; after each X-scan</li> <li>blank line after each Y-scan</li> <li>Some users have gone to I2 format</li> </ul> |
| 9     | ((rho(i),i=1,nx),j=1,ny),k=1,nz)                 | <ul style="list-style-type: none"> <li>Density in each voxel</li> <li>Same format as above</li> </ul>                                                                                                                    |

# Phantom definition: \*.egsphant file


```
4
AIR700ICRU
LUNG700ICRU
ICRUTISSUE700ICRU
ICRPBONE700ICRU
1.00000000 1.00000000 1.00000000 1.00000000
36 45 30
6.98250008 7.48611116 7.98972225 8.49333382 8.99694538 9.50055695
12.0186148 12.5222263 13.0258379 13.5294495 14.0330610 14.5366726
17.0547276 17.5583382 18.0619488 18.5655594 19.0691700 19.5727806
22.0908337 22.5944443 23.0980549 23.6016655 24.1052761 24.6088867
23.4025002 23.9061108 24.4097214 24.9133320 25.4169426 25.9205532
28.4386063 28.9422169 29.4458275 29.9494381 30.4530487 30.9566593
33.4747162 33.9783287 34.4819412 34.9855537 35.4891663 35.9927788
38.5108414 39.0144539 39.5180664 40.0216789 40.5252914 41.0289040
43.5469666 44.0505791 44.5541916 45.0578041 45.5614166 46.0650291
-9.00000000 -8.50000000 -8.00000000 -7.50000000 -7.00000000 -6.50000000 -
-4.00000000 -3.50000000 -3.00000000 -2.50000000 -2.00000000 -1.50000000 -
1.00000000 1.50000000 2.00000000 2.50000000 3.00000000 3.50000000
6.00000000
1111111111113333321111111111111111
1111111111113333332111111111111111
1111111111233333332111111111111111
1111111111233333332211111111111111
1111111111233333332122111111111111
1111111111233433311221111111111111
1111111122123343321222111111111111
1111112211233433432222221111111111
11111233333433433333221111111111
11112333333443333333221111111111
11123333333433333333332211111111
11123333333333333333332111111111
1123333333333333333333443221111111
1123443333333333333333443211111111
1234443333333333333333444321111111
1333444333333333333333443333311111
2333344333333333333333433333211111
2333344433343343343344433333211111
3333344443444333444443443333211111
```

# Isodode contours: dosxyz\_show

- Dose visualization tool based on Motif (Lesstif) libraries
- To use: `dosxyz_show phantom.egsphant [result.3ddose]`







**Joint ICTP-IAEA Workshop on Monte Carlo Radiation Transport  
and Associated Data Needs for Medical Applications**

28 October – 8 November 2024

ICTP, Trieste, Italy

## Lecture 19

---

# DOSXYZnrc calculations with CT input

**Blake Walters**

Metrology Research Centre  
National Research Council Canada

---



Government  
of Canada

Gouvernement  
du Canada

